Microscopio Modelo 701 LED

Cod. 5313111

Manual de Instrucciones / Instructions Manual

INDICE GENERAL

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Content

1. TECHNICAL DATA	2
2. INTRODUCTION, MICROSCOPE STRUCTURE	3
3. UNPACKING.....	3
4. PREPARATION	3
4.1 OBJECTIVES.....	3
4.2 HEAD	3
4.3 CONDENSER.....	3
4.4 FILTER.....	3
4.5 EYE PIECES.....	3
5. OPERATION	3
5.1 AUTOMATIC FOCUSING STOP.....	4
5.2 PUPILLARY DISTANCE ADJUSTMENT	4
5.3 FOCUSING	4
5.4 BLUE DAYLIGHT FILTER.....	4
5.5 MECHANICAL STAGE	5
5.6 OPERATING CONDITIONS.....	5
6. MAINTENANCE	6
a. CLEANING.....	6
b. LAMP REPLACEMENT	6
7. WARRANTY	6


BINOCULAR MICROSCOPE 701 LED

1. TECHNICAL DATA

- Binocular head inclined at 45° and 360° rotation.
- Interpupillary distance adjustable between 52 and 72mm.
- Dioptic adjustment ± 5 dp in both eyepieces, to compensate the differences in visual acuity.
- **Pair of wide field** DIN eyepieces **WF10x** (18Ø)
- **Pair of wide field** DIN eyepieces **WF16x** (11Ø).
- **DIN Achromatic bright field objectives:** **4x** (0.10 NA), **10x** (0.25 NA), **40x** (0.65 NA) retractable and **100x** (1.25 NA) oil immersion and retractable; colour coded for quick identification.
- **Total magnification: 1600x**
- **Quadruple lens nosepiece** with two-way rotation and stop notches.
- **Large stage, 135 x 145mm**, height adjustable and gear mechanism, by bilateral coaxial focussing knobs. **Micrometric** graduated to 0.002mm and **macrometric** with tension adjustment (30mm travel).
- **Displacement limiter** with pre-adjustable safety stop.
- **Mechanical stage** (sliding carriage), with 0,1mm graduation scale, consisting on a smooth rack and pinion mechanism (x-y) **40mm longitudinal and 80mm transversal**.
- **1.25 A.N. Abbe condenser**, with iris diaphragm and filter holder. Height adjustable by rack and pinion.
- **Transmitted illumination (diascopic)**, with removable condenser lens and 1W low-voltage LED lamp with dimmable intensity or fork illumination by means of a flat concave reflecting mirror for optimal illumination by an external source.
- **Very strong stand**, which incorporates in its base the illuminator, intensity potentiometer and power switch.
- **Dimensions:** height 370mm x width 200mm x depth 225mm.
- **Weight:** 14kg.
- **Voltage:** 220V/50Hz
- Supplied with blue, green and yellow filters set, spare fuse, immersion oil bottle, plastic protective cover and wooden box with handle and lock.

2. INTRODUCTION, MICROSCOPE STRUCTURE

The 701 series microscopes are designed for clinical professionals and educators, as well as advanced science, veterinary and medical students.

- | | | | |
|---|-------------------|-------------------------|-----------------------------|
| 1. Eyepieces | 2. Binocular head | 3. Nosepiece | 4. Objectives |
| 5. Slip clutch | 6. Stage | 7. Condenser | 8. Condenser focussing knob |
| 9. Mirror | 10. Base | 11. Body | 12. Focussing lock |
| 13. Coarse focussing knob | | 14. Fine focussing knob | |
| 15. Coarse focussing tension control knob | | | |

3. UNPACKING

Carefully unpack the microscope from the box. Before disposing the packaging, please check the included items against the standard equipment list, to ensure you have received all the components:

1. Base, including support arm, focussing mechanism, objective lens holder and integrated illumination system.
2. Binocular observation head.
3. Achromatic DIN lenses: 4x, 10x, 40x and 100x.
4. Pair of eyepieces: DIN WF 10 x (18 mm $\varnothing$) and WF 16 x (11 mm $\varnothing$).
5. Abbe condenser with iris diaphragm and filter holder.
6. Filters: blue, green and yellow.
7. Dust cover.
8. Oil bottle for immersion.
9. Wooden box with carrying handle and lock.

4. PREPARATION

4.1 OBJECTIVES

To prevent damage to the lens during attachment to the lens holder, turn the coarse adjustment knob counter clockwise for right hand adjustment to lower the lens holder.

4.2 HEAD

Slightly loosen the fixing screen on the head and carefully attach the head to the microscope arm. Lock the head in the desired position by tightening the set screw.

4.3 CONDENSER

Slightly loosen the knurled fixing screw of the condenser and insert the condenser from below. Position the condenser so that the aperture diaphragm knob is conveniently accessible. Tighten the condenser fixing screw with your fingers. To eliminate light irregularities when using low power lenses, such as 4x & 10x, raise and lower the condenser using its adjustment knob.

4.4 FILTER

The filter holder is located in the ABBE condenser.

4.5 EYE PIECES

Place the eyepiece(s) into the eyepiece tubes.

5. OPERATION

Connect the microscope power cord to a properly grounded electrical outlet.

As a precaution, the power cord grounding has been connected to the base of the microscope. Never use an adapter between the power cord and the power supply, as this will disable the described precaution.

Turn ON the illumination switch. To obtain the desired illumination, slide the light controller rheostat to the desired position.

The light intensity degree depends on individual conditions such as the contrast of the specimen, the magnification of the objective, the vision of the individual, etc. It is not recommended to use too much or too little light. The energy-saving illumination indicator lamp remains ON while the illuminator is in use.

5.1 AUTOMATIC FOCUSING STOP

The 701 microscope is equipped with an autofocus locking device to prevent the microscope stage from moving during certain applications. If the autofocus lock is required, focus and then tighten the autofocus knob by turning it clockwise. To ensure free movement of the stage, turn the autofocus lock knob counter clockwise before focussing.

5.2 PUPILLARY DISTANCE ADJUSTMENT

The **701 microscope** is designed to accommodate individual pupillary distances (DEP). To ensure that the DEP is set to a comfortable position, hold the microscope head with both hands and vary the separation between the eyepiece tubes until the entire field is visible simultaneously through both eyepieces.

Turn the dioptre ring on the left eyepiece tube if you are using lenses until the sample image is in sharp focus.

If the interocular and dioptre distances are set correctly:

- the entire field of view can be seen comfortably without moving the head.
- only fine adjustment is required to bring the images into sharp focus by rotating the lenses as they are used.

5.3 FOCUSING

Coarse adjustment –turn counter clockwise the coarse adjustment knob on the right hand. The stage will stop at the bottom of its travel. To raise the slide, turn the knob in the opposite direction.

Place the specimen stage in the centre of the slide. Using a 10x objective and an 10x eyepiece for preliminary observation, focus the specimen by using the coarse adjustment knob.

To prevent damage to the objective lens:

Without using the eyepieces, look directly at the stage and raise the slide by slowly turning clockwise the coarse adjustment knob on the right side. Continue raising the stage until the cover glass is close to touch the objective or until it reaches its upper limit.

While viewing the sample through the eyepieces, turn counter clockwise the coarse adjustment knob to lower the slide and bring the sample into focus.

Adjust the autofocus lock as described in section 5.1. The other objectives shall remain in focal pair without danger of touching the specimen stage. If the slide moves too freely, adjust the tension by turning clockwise the fine adjustment tension knob to tighten or counter clockwise to loosen.

5.4 BLUE DAYLIGHT FILTER

If the light directed to the sample is excessively yellow, place the blue daylight filter in the condenser filter holder.

5.5 MECHANICAL STAGE

When placing a sample on the mechanical slide, the range of travel is approximately 40mm x 80mm. The mechanical slide control knobs, located at the right side and opposite the find adjustment knob, must be turned to change the visual position of the sample on the mechanical slide.

The scale on the upper back section provides two readings:

- The main scale "m" is graduated in millimetres.
- The Vernier scale is represented by "n".


$$1 \text{ "n" division} = 9 \text{ "m" divisions} / 10$$

Note: The picture shows location of "n" between 12mm & 13mm. The closest point to the alignment of "m" and "n" is used to determine the fraction of millimetres to be added.

This picture shows that the 20 of "m" is aligned to the 8 of "n", so it is 0.8mm.

Therefore, the reading is 12mm of "m" plus 0.8mm of "n", that is 12.8mm.

This reading is the result of the "m" indication plus the reading at the point with the closest alignment.

WARNING: Carefully and completely fold the mechanical slide when you are finished using the microscope. If it is fully extended in either directions, the slide may be damaged when attempting to place the microscope in the carrying case.

5.6 OPERATING CONDITIONS

The **701 microscope** can be operated over a wide range of conditions with a minimum adjustment. However, there are variables that must be changed depending on the sample.

The light level should be set at a low point that provides comfortable visibility, thus extending the lamp life considerably. The optimum lamp setting will be change depending on the density of the sample, and the magnification of the objective.

When examining samples with inherently low contrast, you can improve this contrast by reducing the numerical aperture (NA) of the condenser to approximately 75% of the objective's NA. This is best done by removing one of the eyepieces and looking down through the eyepiece tube into the illuminated rear aperture of the objective. While looking at this aperture, close the condenser diaphragm aperture until the size of the rear aperture of the lens is reduced to approximately 75% of its original size.

Samples whose contrast is so low as to be invisible, or barely visible, can often be seen more clearly by closing the aperture diaphragm of the condenser to the size of a pinhole. An alternative and often superior method of improving contrast is to use oblique illumination. This is achieved by covering 1/2 or more of the condenser aperture with opaque paper.

Both methods of contrast enhancement introduce diffraction effects that reduce the microscope's ability to highlight fine details of the sample and it is not recommended in cases where critical examination of the sample is required. In that case, use the phase contrast microscopy technique.

6. MAINTENANCE

Keep the equipment clean and with a minimal routine maintenance to ensure a long life for the microscope.

Important: Always protect the microscope with a dust cover when not in use. Always keep the objectives and eyepieces in place when storing the microscope.

a. CLEANING

The external surface should be cleaned from time to time with a soft damp cloth. Only use mild soap solutions to clean more heavily damaged areas.

External surfaces of the optics should be inspected for dust and dirt. Use a compressed air bulb to blow dust off the optics. If necessary, use a soft cloth or cotton toothpick, moistened with water and a mild soap solution to clean the optics surfaces. As a last resource, use alcohol or xylene.

Note: Avoid using solvents, as they can cause problems with adhesive with which the optics are attached. Also, solvents can attract grease from the frames, making cleaning more difficult.

Each time immersion oil is used, wipe the objective and the slide with a soft cloth or lens tissue paper. Occasionally, the underside of the mechanical slide may drip with oil and require cleaning.

Disassemble the mechanical slide by removing the three screws on the underside of the slide table. When wiping the oil, be careful not to remove the slide rails or contaminate them with lubricant. Reinstall the mechanical slide carefully to avoid damage to the slide rails, rack and pinion.

b. LAMP REPLACEMENT

The only component that may require periodic replacement in the microscope is the 1W LED lamp.

Note: The lamp glass surface may be extremely hot. Do not attempt to change the lamp before it has been allowed to cool completely or without wearing a suitable skin protection.

Before replacing the lamp, switch the equipment off. Place the mechanical slide in the extreme right position and carefully place the microscope on its left side. Open the lamp door on the bottom plate by pulling the knob.

After the lamp has cooled down, carefully remove it from the lamp holder. Put the new lamp in place and close the door. Never operate the microscope illuminator unless the lamp door is properly closed.

7. WARRANTY

This product has a warranty period of one year. The warranty does not cover lamps or damages caused by an improper use or causes beyond COMECTA, S.A.'s control.