REFRACTÓMETRO CLÍNICO 2 ESCALAS

Modelo C-6 Cód.5901006

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Model: C-6 Part number: 5901006

OPERATING MANUAL

Read the instructions carefully before using the equipment.


Figure 1

DESCRIPTION

Precision optical instrument for medical applications.
Stands out for its small volume, light weight and to be easily portable.

APPLICATIONS

1. Albumin in serum measurement.
2. Urine specific density

TECHNICAL DATA

- | | |
|-------------------------|-----------------------|
| 1. Magnification: | 3.4X |
| 2. Measurement range: | |
| Albumin in serum: | from 0 to 12 g/dl |
| Graduation: | 0,2 g/dl |
| Urine specific density: | from 1,000 to 1,040sg |
| Graduation: | 0,002 |
| 3. Measurements: | 27 x 40 x 160 mm. |
| 4. Weight: | 176 g. |

OPERATION

1. Point the refractometer towards the light and adjust the dioptic ring (5) until graduation is clearly read.
2. Adjust the balance point 0:
Open the cover (2), and place one or two drops of distilled water on the prism surface, close the lid, and adjust (3) in order that the light-dark contrast coincides with W line.
3. Open the cover again (2) and clean the prism and the cover surface with a soft cloth. Place 1 or 2 drops of serum (plasma), close the lid with a slight pressure. Read the graduation which, with light / dark position, will directly indicate the concentration albumin in serum.
4. When measuring the urine specific density, you must proceed in the same way.
5. After measurements, always clean the cover and the prism.

MAINTENANCE

1. The adjustment with distilled water and the sample reading should always be done at the same temperature. If this changes significantly, the balance point should be adjusted once every 30 minutes.
2. After its use, do not use water to clean the instrument.
3. As it is a precision optical instrument, it should be handled with special care and attention. Do not touch the optical surfaces and keep it in a dry place.
4. If you follow these use conditions, the duration of the instrument is guaranteed. The optical elements cannot be changed.
5. The specific weight read must be increased in 0.005 units when the added protein is of about 1% and 0.0004 units if the added sugar is 1%.