
Instructions Operating

Read the instructions completely before operating this set.

1. Prism 2. Cover plate 3. Correct screw 4. Mirror tube 5. Eyepiece (Adjusting ring of diopter)

Figure 1

APPLICATION

- * These models are applicable to test the concentration of sugar in solution
- * Models C4 and C5 and the ATC versions are with a higher precision for testing the sugar solution at low concentration, such as: fruit juice, soft drink, grape wine etc.
- * Model C2 and the ATC version is applicable for general purpose of sugar testing.
- * Model C3 and the ATC version are applicable especially to test sugar concentration in fruit juice, condensed milk, liquid state sugar and jam, etc.

SERIES


STYLE	MODEL	MEASURING RANGE	ACCURACY	SIZE DIAM X L/mm	WEIGHT
BRIX	C4	0 - 10%	0,1%	29 x 160	175 g
	C5	0 - 20%	0,1%	29 x 160	175 g
	C2	0 - 32%	0,2%	29 x 160	175 g
	C3	28 - 62%	0,2%	29 x 148	175 g

OPERATION PROCESS

1. Turn the inclining prism of refractometer upside to the direct light source and put one eye close to the eyepiece, and then adjust the diopter by turning eyepiece till the scale with digits and letters inside emerges from the image clearly.
2. Calibration before using: lift the cover plate (2) and drop one or two drops of pur distilled water on the surface of the prism, then close the cover plate on the surface of prism tightly. Adjust the calibration screw (3) till the light/dark boundary meets the zero line.
 - Model C3 is calibrated by a saturated sodium chloride solution, The refractive index of he saturated sodium chloride solution is: adjust at 29.9% for 15°, adjust at 29,5% for 20°, adjust at 29,2% for 25°
3. Open the cover plate (2) and clean the surface of prism and cover plate by soft cotton flannel, then drop 1-2 drops of the liquid needs to be tested on the primsm surface. Close the cover plate to the prism surface loghtly, the corresponding testing value can be read from the light/dark boundary crossed on scale.
4. After cleaning and drying the surface of the prism and cover plate by moist gauze, store the refractometer at an appropriate place.

ATTENTION & MAINTENANCE

1. The calibration and testing should be at the same temperature. If the temperature is varied, the refractometer should be re-caibrated per each 30 minutes.
2. The prism should be cleaned completely in case that the residual impurity would cause error in testing result. Don't use water to wash refractometer in order to avoid the built-in components being soaked and damaged.
3. Storage: keep it at a dry, clean and ventilated place and keep it away from the corrosive & volatile material and dampness in case of mildewing and fogging on optical glass.
4. Avoiding touchig, scratching on the optical surfaces and avoiding strong shock during transportation. Handling softly and gently.
5. The refractometer should be operated in accordance with this instruction manual strictly so as to keep in good condition and good optical performance.

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