MODELO 'NDJ-1'

COD. 5120230

VISCOSÍMETRO ROTACIONAL
ROTATING VISCOMETERS

INFORMACIÓN GENERAL**GENERAL INFORMATION**

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PACKING LIST

The standard equipment consist of the following components:

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APPLICATIONS

The rotational viscosimeter model 'NDJ-1' is a novel instrument to be used in measuring the viscous resistance and the absolute viscosity of fluid.

It is widely applicable to determine the viscosity of various kinds of fluids, such as grease, paint, plastics, food products, pharmaceuticals, bonding agent, etc.

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MAIN TECHNICAL DATA

Range:	Medium viscosity 10-100.000mPa.s Medium viscosity 4 rotors of Nos. 1 to 4.
Speed of rotor:	6rpm, 12rpm, 30rpm and 60rpm
Measur.tolerance:	+/-5% (Newtonian fluid)
Power supply:	220V +/-10V.50Hz
Dimensions:	300 x 300 x 450mm
New weight:	1.5kg (excluding support stand)

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PRINCIPLE OF STRUCTURE


1. As shown in Fig.1, the synchronous, which is connected to the graduated disc, rotates in a steady speed, and drives the rotor through the hairspring and the shaft. If there is no resistance on the fluid, the rotor will rotate synchronously with the hairspring, indicator and graduated disc, the indicator thus indicated will be at the reading "0". On the contrary, if the rotor is subjected to the viscous resistance of fluid, there will be a torque at the hairspring. After it is balanced against the viscous resistance, the indicator connected to the hairspring will then indicate a particular reading (the torsional moment of the hairspring) on the graduated disc. Multiplying of this reading with the specific coefficient will obtain the viscosity of the fluid (mPa.s).

2. According to the requirements in determining the viscosity of the fluid, the speed of the viscosimeter can be adjusted into four speeds by a knob.

3. Four rotors with a protector rotor device will be provided.

The user can select the appropriate rotor coupled with the appropriate speed according to the viscosity of the fluid to be measured.


4. The viscosimeter is provided with a indicator fixing control device, so that when the speed is too fast (30rpm, 60rpm), that reading is impossible during rotation, the control device can be pressed to stop the indicator for easy reading.

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
5. The protection rack is used to stabilize the measurement and to protect the rotor.
6. The viscosimeter can be operated hand hold. As it is provided with the support stand and elevation device, it is recommended to fix the viscosimeter on the support stand when making small amount and fixed temperature measurement.

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INSTALLATION

1. Take out the packing box, support stand, column, etc. from the case.
2. Screw the column into the hole at the rear of the support stand with the teethed side facing the front of the support stand. Tighten the nut with wrench to prevent the column rotating.
3. Turn the knob of the elevation clamp, and check the flexibility and selflocking of the elevation clamp. If it is found to be too loose or too tight, adjust the stop screw on the clamp with a wrench for socket head cap screw. It should be moved up and down freely, but a little tight is preferable to prevent the viscosimeter dropping down after it is installed.
4. Open the packing box, take out the rotor stand and the viscosimeter, install it on the elevation clamp, tighten the fixing screw on lever, connect it to the power supply, remove the rubber band on the indicator control lever, release the yellow screw at the lower end of the instrument and lastly remove the yellow packing ring.
5. Regulate the levelling screw to keep the instrument level.

OPERATION

1. Prepare the fluid to be measured, put in a beaker or cylindrical container with its diameter not less than 70mm and control the temperature of the fluid precisely.
2. Install the protection rack (Turn right to install and left to remove).
3. Screw the selected rotor into the connecting worm (Turn left to install and right to remove). Turn the elevation knob, slowly lower the instrument, the rotor will gradually dip into the fluid to be measured until the fluid surface mark on the rotor is level with the fluid surface, then adjust the levelling of the instrument. Press indicator control rod, switch on the motor, turn the speed changing knob making the required speed upward to align with the indication point for speed, release the indicator control rod causing the rotor to rotate in the fluid.

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After several rotation (generally 20-30 seconds) until the indicator tends to be stable, press indicator control rod (Caution: 1. Do not use excessive force. 2. When the speed is quite slow, it is not necessary to use control rod. The readings can be directly made), causing the reading to settle, switch off the motor, the indicator will stay at the readout window and then reading can be taken. After the motor has switched off but the indicator does not stay at the readout window, press the indicator control rod continuously and repeatedly switch on and switch off the motor. After several practices the user will skillfully manage the technique of keeping the indicator to stay in the window for reading.

4. When the value indicated by the indicator is too high or too low, change the rotor or speed so the reading is within the preferable range of 30-90 graduation.

5. Selection of range, coefficient and rotor speed:
 - A. First estimate the approximate range of viscosity of the fluid to be measured, then select the appropriate rotor and speed according to Table of Range.

For instance if the fluid to be measured is approximately 3.000 mPa.s the following combinations can be selected:

Rotor no.2 6 rpm
 or Rotor no.3 30rpm

- B. When the approximate viscosity of the fluid to be measured is difficult to estimate, suppose a higher viscosity first, then try the rotor of lower number to higher number and lower speed to higher speed. The principle is that for the fluid of high viscosity select small rotor and slow speed, for the fluid of low viscosity select large rotor and fast speed.

- C. Coefficient: The reading indicated by the indicator on the graduated disc should multiply by the particular coefficient in Table or Range to get the absolute viscosity (mPa.s), i.e.

$$\eta = k \cdot a$$

where

η is absolute viscosity

k is coefficient

a is the reading indicated by the indicator (deflected angle)

- D. Correction of frequency error: When the frequency of power supply is not correct, it can be corrected according to the following formula.

$$\text{Practical visc.} = \text{Indicated visc.} \times \frac{\text{Nominal freq.}}{\text{Practical freq.}}$$

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POINTS FOR ATTENTION

1. *This instrument should be used in room temperature.*
2. *The measurement should be made within the specified frequency and voltage, or the precision of the measurement will be affected.*
3. *Support the instrument on the support stand preferably. If it is operated hand held, keep the instrument steady and level.*
4. *Care should be taken in installing and removing the rotor. To install and remove the rotor the connecting worm should be slightly lifted. Do not use excessive force and apply lateral force on the rotor to prevent flexure of the rotor.*
5. *When the rotor is installed, the instrument should not be placed on its side or upside down.*
6. *Before pressing the indicator control rod, the motor should not be switched on. When the motor works, the speed can be changed.*
7. *The interface of the connecting worm and the rotor as well as the thread should be kept clean, or it will affect the correct connection of the rotor and the stability during rotation.*
8. *When lifting or lowering the instrument, hold it with hand to prevent dropping due to its dead weight.*
9. *After each operation, clean the rotor immediately (do not clean the rotor when it is on the instrument). After cleaning, put the rotor on the rotor rack properly.*
10. *Do not dismantle or adjust the parts of the instrument arbitrarily.*
11. *When moving or transporting the instrument, use a rubber band to wrap the indicator control rod, sleeve it into the yellow packing ring and lift the connecting worm, then tighten with screw.*
12. *Most of the high viscosity fluid, such as suspension, emulsion, high polymer, etc. are usually "Non-Newtonian fluid", their variation of viscosity is according to the shear rate and change of time. Therefore measurement under different rotor, speed and time may give nonuniform results. This is normal, and it is not the fault of the instrument. (Usually when measuring non-Newtonian fluid, rotor, speed and time should be specified).*

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13. To obtain a precise value of measurement, the following points should be observed:

- A-** Precisely control the temperature of the fluid to be measured.
- B-** Dip the rotor into the fluid to be measured for a sufficient time to warm up, so that it will have same temperature as the fluid.
- C-** Ensure the uniformity of the fluid.
- D-** Put the rotor in the center of the container during the measurement.
- E-** Avoid bubbles adhering under the rotor when the rotor is dipped into the fluid.
- F-** When changing the rotor and the speed, keep the reading on the graduated disc a little higher.
- G-** Use the protection rack to measure.
- H-** Ensure the cleanliness of the rotor.
- I** When the frequency of the power supply is not correct, use the above mentioned formula for correction.
- J** Strictly adhere to the operation instructions for measurement.

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GUARANTEE

This product is guaranteed for one year. The guarantee does not cover damage caused by incorrect use or causes beyond the control of COMECTA, S.A.

Any manipulation of the apparatus by personnel not authorized by COMECTA, S.A. automatically cancels the guarantee.