

Automatic Module

Hei-VOLUME Distimatic

The clever way of
continuous evaporation

With the unique residue drainage
of the Distimatic automatic module,
unlimited quantities can be automatically
distilled in unattended, continuous operation.

**24 hours a day,
7 days a week.**



- Automatic supply and removal of solvent mixture and condensate
- Almost one hundred percent solvent separation in automated unattended and continuous operation
- Individual determination of the concentration of the final product in order to prevent drying-out, crystallization and powder formation
- Automated supply of small volumes to avoid foaming
- Exchange of the storage and collecting vessels possible during the running process
- Permanent process monitoring via sensors, incl. full system shutdown in case of an emergency
- Compatible with all Heidolph rotary evaporators and large-scale rotary evaporators

Hei-VOLUME Distimatic

Automatic 24/7 Evaporation

Functional principle

- 1 Vacuum valve
- 2 Collector valve
- 3 Ventilation valve
- 4 Condensate pump
- 5 Residue pump



Distimatic Integration Package

For integration into an existing vacuum and cooling liquid system without drainage of residue

- Hei-VAP Ultimate Control ML G9B XL
- Tubing set
- Distimatic without automatic drainage of residue

P/N 576-82915-00

Distimatic Integration Package 24/7

For integration into an existing vacuum and cooling liquid system including drainage of residue

- Hei-VAP Ultimate Control ML G9B XL
- Tubing set
- Distimatic incl. automatic drainage of residue

P/N 576-83915-00

Distimatic Workstation

Complete package without drainage of residue

- Hei-VAP Ultimate Control ML G9B XL
- Hei-CHILL 1200
- Kryo 30 (20 l)
- RS 232 cable
- Rotavac Vario Control
- Condenser for Rotavac Vario Control
- Tubing set
- Distimatic without automatic drainage of residue
- Heating bath liquid (5 l)
- 5 l Evaporating flask

P/N 576-80915-00

Distimatic Workstation 24/7

Complete package including drainage of residue

- Hei-VAP Ultimate Control ML G9B XL
- Hei-CHILL 1200
- Kryo 30 (20 l)
- RS 232 cable
- Rotavac Vario Control
- Condenser for Rotavac Vario Control
- Tubing set
- Distimatic incl. automatic drainage of residue
- Heating bath liquid (5 l)
- 5 l Evaporating flask

P/N 576-81915-00

Discover more about our new Automatic Modul - Distimatic under www.heidolph-instruments.com